

"YOUR MOBILITY MADE EASY!"

"OT NEWS YOU CAN USE" ...

AEM TIP: Always remember when your mobility scooter or powered wheelchair is not in use to keep your batteries fully charged. As flat uncharged batteries left for a period of time may not take a charge and will need to be replaced. Your charger will keep your batteries fully charged when not in use and it wont over charge them. So when not in use for a period of times always have your mobility scooter or powered wheelchair charged regularly i.e. once a week.

NEW PRODUCT LAUNCH

The Rascal Partner

This deluxe Community Mobility Scooter is a robust, manoeuvrable mobility scooter for convenient indoor and outdoor use. The Partner mobility scooter is of fixed frame construction and not designed to be fully dismantled for transportation, however it is compact enough to fit into some larger cars, jeeps and vans. It is ideal for use around your community and with its larger wheels and seat your comfort is guaranteed. There is two models available the standard which will travel 8 miles and the deluxe which will travel up to 16 miles on one fully charged battery.



Specifications

Range: 15 miles	Max User Weight: 17.5 stone or 113kg
Length: 111.8 cm (44")	Width: 60cm (24")Height: 91.7cm (36.1")
Weight of Batteries: 43.9kg or 6.9 stone	Ground Clearance: 5cm (2")
Max Safe Slope: 8 degrees	Max speed: 4mph
Front wheel size: 20cm x 6.7cm	Rear wheel size: 23.8cm x 7cm
Turnaround width: 1.8m or 5.9ft	Full 12 month warranty

COSYFEET EXTRA ROOMY SHOES, SLIPPERS & SANDALS NOW IN STOCK

Cosyfeet are in stock now, contact us to arrange a fitting of these shoes, sandals and slippers in the comfort of your own home. These are ideal for people who suffer with swollen or bandaged feet and are light weight. We stock the most popular sizes and styles, if you require something we do not have in stock normal order time is just a few days.



Sale now on . If you have any questions or if you are finding it hard to source any products or services please feel free to call us and we will be more than happy to help you in any we can.

**CONTACT US NOW
014553168**

ADVANCE ELECTRICAL MOBILITY

Sales (Hire*)

Maintenance

Repair

FULL RANGE &
CATALOGUES OF

MOBILITY
SCOOTERS*
POWERED WHEEL-
CHAIRS*
RISER RECLINER
ARMCHAIRS
PROFILING BEDS
BATHLIFT SEAT
PATIENT LIFTING
HOIST*
STAIRLIFTS
MANUAL WHEEL-
CHAIRS*
WALKING AIDS*
HIGH BACK
ARMCHAIRS

ADVANCE ELECTRICAL
MOBILITY

Sales (Hire*)

Maintenance

Repair

4 CRUMLIN VILLAGE
DUBLIN 12

Phone: 01-4553168

Fax: 01-4550763

E-mail:

info@aemobility.com

Web:

www.aemobility.com

